


## DIAZOAMINO BENZENE



### Reagents:

conc. HCl – 6.7 ml  
Aniline – 4.6 ml  
 $\text{NaNO}_2$  – 1.75 g  
Sodium acetate – 7 g

### Apparatus:

Erlenmeyer flask with ground glass 200 ml  
Dropper 25 ml  
High beaker 50 ml  
glass funnel  
thermometer up to  $150^\circ\text{C}$   
crystalliser  
glass rod  
Petri glass

### Preparation:

In a Erlenmeyer flask, 25ml of water, 6.7ml of concentrated hydrochloric acid and 4.6ml of aniline are placed. The flask is shaken vigorously and 17g of broken ice is added.

A 30% aqueous solution of sodium nitrate (III) (5ml) is then added over 5-10 minutes, shaking continuously. The reaction mixture is left to stand for 15 minutes (shaking quite frequently) and then a 35% aqueous solution of sodium acetate (21ml) is added within 5 minutes.

The yellow diazoaminobenzene precipitate is shaken at below  $20^\circ\text{C}$  (ice is added if necessary) over 45 minutes, then drained, washed with about 80ml of cold water (the precipitate should be pressed as thoroughly as possible) and air-dried

Melting point:  $90-91^\circ\text{C}$